

IRON (CAB)

INTENDED USE

For the quantitative determination of iron in human serum or heparinized plasma.

CLINICAL SIGNIFICANCE

The iron is the component of a great number of enzymes. The myoglobin, muscular protein, contains iron, as well as the liver. Iron is necessary for the hemoglobin production, molecule that transports oxygen inside red globules. Their deficit in the last causes the ferropenic anemia. High levels of iron are found in hemochromatosis, cirrhosis, hepatitis and in increased transferrin levels. The variation day to day is quite marked in healthy people. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

METHOD AND PRINCIPLE

Serum iron reacts with chromazurol-B and CTMA-Br to form a complex deeply blue colored. The absorbance measured at 630 nm is directly proportional to the amount of iron in the sample.

The intensity of the color formed is proportional to the iron concentration in the sample.

REAGENT COMPOSITION

R 1 Acetate Buffer pH 4.9 100 mmol/L

R 2 Color CAB 40 mmol/L

IRON CAL Iron aqueous primary standard 200 µg/dL

WARNINGS AND PRECAUTIONS

1. Iron STD: proceed carefully with this product because due its nature it can get contaminated easily.
2. It is recommended to use disposable material. If glassware is used the material should be soaking for 6 h in diluted HCl (20% v/v) and then thoroughly rinsed with distilled water and dried before use.
3. Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
4. Use clean disposable pipette tips for its dispensation.
5. The reference values are strongly method dependent.

STORAGE AND STABILITY OF REAGENT.

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

REAGENT DETERIORATION

Presence of particles and turbidity.

Blank absorbance (A) at 630 nm > 0.500 abs.

SPECIMEN COLLECTION AND STABILITY

Serum or heparinized plasma. Free of hemolysis and separated from cells as rapidly as possible. Stability of the sample: 2-8°C for 7 days.

INTERFERENCES

Hemolyzed samples are rejected, since erythrocytes contain iron and therefore falsely elevate the serum results. A list of drugs and other interfering substances with iron

determination has been reported.

ASSAY PROCEDURE FOR SEMIAUTO ANALYZER.

Wavelength : 630 nm

Temperature :37 deg c

	Blank (µl)	Std (µl)	Test (µl)
Reagent	1000	1000	1000
Std	--	50	--
Test	--	--	50

Mix and incubate the tubes for 10 minutes at 37 deg c. Measure the optical density of standard against reagent blank to generate factor . Multiply test absorbance with the factor to get the IRON results as below.

CALCULATIONS

$(A) \text{ Sample} - (A) \text{ Reagent Blank} \times 200 \mu\text{g/dL}$

$(A) \text{ Standard} - (A) \text{ Reagent blank}$

QUALITY CONTROL

Use control sera with known normal and abnormal values to monitor the integrity of the reaction. Values should be those acceptable for this method and temperature.

EXPECTED VALUES

Male 65 - 175 µg/dL

Female 40 - 150 µg/dL

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 0.850 µg/dL to linearity limit of 500 µg/dL. If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision

	Intra-assay (n=20)		Inter-assay (n=20)	
Mean	113	250	111	249
SD	0.89	0.72	3.51	6.29
CV	0.79	0.29	3.17	2.52

Sensitivity: 1µg/dl

Results obtained using competitors reagents did not show systematic differences when compared with other commercial reagents. The results obtained using 50 samples were the following : Correlation coefficient (r)2:0.9934.

Regression equation : $y = 1.0243x - 3.877$. The results of the performance characteristics depend on the analyzer used.

GENERAL TECHNICAL PARAMETERS

Mode	End point.
Wavelength	630 nm
Reaction direction	Increasing
Reagent Blank	Yes
Sample Vol	0.050 ml
Reagent Vol	1.000 ml
Measuring time	10 min at 37 deg c
Reagent Blank Abs Max	NMT 0.500 Abs
Calibration Mode	1-point
Standard Conc	200 µg/dL
Linearity	500 µg/dL

Decimal Places	1
Unit	µg/dL
Ref Low	Male 65 µg/dL /Female 40 µg/dL
Ref high	Male 175 µg/dL /Female 150 µg/dL

REFERENCES

1. Perrotta G. Iron and iron-binding capacity. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1063-1065.
2. Itano M M D. Cap Serum Iron Survey 1978 (70): 516-522.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.