


**Rapid test for
FSH**




**For Qualitative
Detection of
FSH (Follicle
Stimulating
Hormone)
Urine Rapid Test**

- Urine samples may be used for testing
- Results available in only 3 minutes
- Easy to use
- Results can be read visually
- No need for an analyser
- Semi-quantitative point-of-care Test for follicle stimulating hormone (FSH)




READING TEST RESULTS



Add 3 drops of (120µL) urine in sample well "S" using the dropper provided. Avoid trapping air bubbles in the specimen well.




Result



Positive

Negative



Invalid