

First View फर्स्ट व्यू


Rapid test for H. Pylori-Ab

For Qualitative Detection of H. Pylori-Ab Whole Blood Rapid Test

- Rapid test for the qualitative detection of IgG antibodies specific to Helicobacter pylori in human whole blood and serum.
- Useful aid in the diagnosis of H. pylori infection
- Easy to use
- CLIA Waived for whole blood
- High accuracy and reliability



READING TEST RESULTS



Add 20µl of serum/plasma/whole blood to the specimen well (S) of the test cassette.



Then add 2 drop of buffer and start the timer.



Result

