

Firstview[®]

One Step Rapid Test Device

Bioline
DIAGNOSTICS
Your safety. Our commitment.

Rapid test for Syphilis


**For Qualitative
Detection of
Syphilis
Whole Blood
Rapid Test**

- Syphilis is a sexually transmitted infection (STI) caused by the spirochete bacterium *Treponema pallidum*.
- Run with PLASMA or SERUM
- Rapid test result in just 15 minutes
- Semi-quantitative point-of-care Test for Syphilis/VDRL

READING TEST RESULTS




Add 1 drops (25 μ L) of serum or plasma in sample well "S" using the dropper provided.




Add 1 drop of the Assay buffer



Result



Positive



Negative



Invalid



ISO 9001
CERTIFIED



ISO 13485
CERTIFIED

Bioline
DIAGNOSTICS
Your safety. Our commitment.

Marketed By

Bioline Diagnostics LLP

H-1478, DSIDC, Narela Industrial Area, Delhi - 110 040

T. +91 11 4170 8100 | E. info@biolinediagnostics.com

www.biolinediagnostics.com